

## GENERAL INFORMATION

**Product name :** Conical head clip, plastic

**Description :** Manufactured in POM (Polyoxymethylene). Will not scratch glass surfaces

## TECHNICAL DATA

reference	type	pcs/pack
CHCP-014-010	orange, for 14/23 joints	10
CHCP-019-010	blue, for 19/26 joints	10
CHCP-024-010	green, for 24/29 joints	10
CHCP-029-010	red, for 29/32 joints	10

## PACKAGING AND LOGISTICS

reference	vol (l)	kg	TARIC	GTIN
CHCP-014-010	0,442	0,025	39269097	08434868013759
CHCP-019-010	0,3	0,036	39269097	08434868013766
CHCP-024-010	0,384	0,052	39269097	08434868013773
CHCP-029-010	1,188	0,076	39269097	08434868013797

## PRODUCT PHOTO



## INFORMACIÓN GENERAL

**Nombre del producto :** Pinzas de plástico para esmerilados**Descripción :** De POM (Polioximetileno). No rayan el vidrio

## DATOS TÉCNICOS

referencia	tipo	unidades por ref.
CHCP-014-010	naranja, para bocas 14/23	10
CHCP-019-010	azul, para bocas 19/26	10
CHCP-024-010	verde, para bocas 24/29	10
CHCP-029-010	rojo, para bocas 29/32	10

## EMBALAJE Y DATOS LOGÍSTICOS

referencia	vol (l)	kg	TARIC	GTIN
CHCP-014-010	0,442	0,025	39269097	08434868013759
CHCP-019-010	0,3	0,036	39269097	08434868013766
CHCP-024-010	0,384	0,052	39269097	08434868013773
CHCP-029-010	1,188	0,076	39269097	08434868013797

## FOTO DEL PRODUCTO



## INFORMATIONS GÉNÉRALES

**Nom produit : Pince de rodage en plastique****Description :** Fabriquée en POM (Polyoxyméthylène). Ne raye pas le verre

## DONNÉES TECHNIQUES

référence	type	unités par ref.
CHCP-014-010	orange, pour col 14/23	10
CHCP-019-010	bleu, pour col 19/26	10
CHCP-024-010	vert, pour col 24/29	10
CHCP-029-010	rouge, pour col 29/32	10

## EMBALLAGE ET LOGISTIQUE

référence	vol (l)	kg	TARIC	GTIN
CHCP-014-010	0,442	0,025	39269097	08434868013759
CHCP-019-010	0,3	0,036	39269097	08434868013766
CHCP-024-010	0,384	0,052	39269097	08434868013773
CHCP-029-010	1,188	0,076	39269097	08434868013797

## PHOTO PRODUIT



## INFORMAZIONE GENERALE

**Nome del prodotto :** Clips in plastica per giunti conici smerigliati

**Descrizione :** In POM (poliossimetilene). Non raschiano il vetro

## DATI TECNICI

referenza	tipo	unità per ref.
CHCP-014-010	arancione, per coni 14/23	10
CHCP-019-010	blu, per coni 19/26	10
CHCP-024-010	verde, per coni 24/29	10
CHCP-029-010	rosso, per coni 29/32	10

## IMBALLAGGIO E DATI LOGISTICI

referenza	vol (l)	kg	TARIC	GTIN
CHCP-014-010	0,442	0,025	39269097	08434868013759
CHCP-019-010	0,3	0,036	39269097	08434868013766
CHCP-024-010	0,384	0,052	39269097	08434868013773
CHCP-029-010	1,188	0,076	39269097	08434868013797

## FOTO DEL PRODOTTO

