

## RECOMMENDED TIPS FOR EASY MICROPIPETTES

	Capacity	Standard Line	Premium Line	Premium Line ULR	Filtered Premium Line
<b>EASY 30+</b>					
EASC-10P-001	100 µl	TIPP-200-***	TIAD-200-***	TIAU-200-***	TIAF-100-***, TIAF-200-***
EASC-20P-001	200 µl	TIPP-200-***, TIPP-201-***	TIAD-200-***, TIAD-300-***	TIAU-200-***	TIAF-200-***
EASC-50P-001	500 µl	TIPP-1K0-***, TIPP-1K1-***	TIAD-1K0-***	TIAU-1K0-***	TIAF-1K0-***
EASC-1KP-001	1000 µl	TIPP-1K0-***, TIPP-1K1-***	TIAD-1K0-***	TIAU-1K0-***	TIAF-1K0-***
<b>EASY 40+</b>					
EASD-P02-001	0,1 - 2,5 µl	TIPP-010-***, TIPP-011-***	TIAD-010-***	TIAU-010-***	TIAF-010-***
EASD-P10-001	0,5 - 10 µl	TIPP-010-***, TIPP-011-***	TIAD-010-***	TIAU-010-***	TIAF-010-***
EASD-P20-001	2 - 20 µl		TIAD-020-***	TIAU-020-***	TIAF-020-***
EASD-P50-001	5 - 50 µl	TIPP-200-***, TIPP-201-***	TIAD-200-***	TIAU-200-***	TIAF-100-***
EASD-10P-001	10 - 100 µl	TIPP-200-***	TIAD-200-***	TIAU-200-***	TIAF-100-***, TIAF-200-***
EASD-20P-001	20 - 200 µl	TIPP-200-***, TIPP-201-***	TIAD-200-***, TIAD-300-***	TIAU-200-***	TIAF-200-***
EASD-1KP-001	100 - 1000 µl	TIPP-1K0-***, TIPP-1K1-***, TIPP-1K2-***	TIAD-1K0-***	TIAU-1K0-***	TIAF-1K0-***
EASD-5KP-001	500 - 5000 µl		TIAD-5K0-***		
EASD-10K-001	1000 - 10000 µl		TIAD-10K-***		
<b>EASY 40 ELITE</b>					
EASD-E02-001	0,1 - 2,5 µl	TIPP-010-***, TIPP-011-***	TIAD-010-***	TIAU-010-***	TIAF-010-***
EASD-E10-001	0,5 - 10 µl	TIPP-010-***, TIPP-011-***	TIAD-010-***	TIAU-010-***	TIAF-010-***
EASD-E20-001	2 - 20 µl	TIPP-200-***, TIPP-201-***	TIAD-200-***	TIAU-200-***	TIAF-100-***
EASD-E50-001	5 - 50 µl	TIPP-200-***, TIPP-201-***	TIAD-200-***	TIAU-200-***	TIAF-100-***
EASD-10E-001	10 - 100 µl	TIPP-200-***, TIPP-201-***	TIAD-200-***, TIAD-300-***	TIAU-200-***	TIAF-200-***
EASD-20E-001	20 - 200 µl	TIPP-200-***, TIPP-201-***	TIAD-200-***, TIAD-300-***	TIAU-200-***	TIAF-200-***
EASD-1KE-001	100 - 1000 µl	TIPP-1K0-***, TIPP-1K1-***, TIPP-1K2-***	TIAD-1K0-***	TIAU-1K0-***	TIAF-1K0-***
EASD-5KE-001	1000 - 5000 µl		TIAD-5K0-***		
<b>EASY W</b>					
EASW-010-001	0,5 - 10 µl	TIPP-010-***, TIPP-011-***	TIAD-010-***	TIAU-010-***	TIAF-010-***
EASW-050-001	5 - 50 µl	TIPP-200-***, TIPP-201-***	TIAD-200-***	TIAU-200-***	TIAF-100-***
EASW-300-001	30 - 300 µl	TIPP-200-***, TIPP-201-***	TIAD-300-***, TIAD-200-***	TIAU-200-***	TIAF-200-***
EASW-1K0-001	100 - 1000 µl	TIPP-1K0-***, TIPP-1K1-***, TIPP-1K2-***	TIAD-1K0-***	TIAU-1K0-***	TIAF-1K0-***